# firesafe™ nozzle

instructions for use

---

5	<b>EN</b>	Firesafe™ Nozzle - Instructions for Use
6	<b>DE</b>	Firesafe™ Nozzle - Bedienungsanleitung
7	<b>FR</b>	Firesafe™ Nozzle - Mode d'emploi
8	<b>ES</b>	Firesafe™ Nozzle - Instrucciones de uso
9	<b>PT</b>	Firesafe™ Nozzle - Instruções de utilização
10	<b>IT</b>	Firesafe™ Nozzle - Istruzioni per l'uso
11	<b>NL</b>	Firesafe™ Nozzle - Gebruiksaanwijzing
12	<b>EL</b>	Firesafe™ Nozzle - Οδηγίες Χρήσης

---

EN Figure 1

DE Abbildung 1

FR Figure 1



ES Figura 1

PT Figura 1

IT Figura 1

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## Intended Use

The Firesafe Nozzle™ is a thermal fuse designed to extinguish an oxygen delivery tube fire by stopping the flow of oxygen if the tube is accidentally ignited and the fire burns back and reaches the Firesafe™ Nozzle.

Oxygen fires are relatively common, often caused by patients smoking whilst receiving oxygen therapy. Oxygen is not flammable but its presence will drastically increase the rate and severity of combustion. Things that will not burn in air may do so in an oxygen enriched environment. Should a fire start in the oxygen delivery tube it will burn back towards the oxygen source, typically an oxygen cylinder, hospital pipeline with a flowmeter or an oxygen concentrator or liquid supply.

It is recommended that two Firesafe devices are fitted to home oxygen delivery circuits. The first should be positioned close to the source of supply and the second close to the patient as shown in figure 1. It is recommended to protect all oxygen supply sources, including hospital oxygen flowmeters and oxygen cylinders, with at least one Firesafe device.

The Firesafe™ Nozzle is a risk control measure intended to be permanently fitted to medical devices that deliver oxygen to patients requiring supplemental oxygen, whether in the home, healthcare facility or in the care of the emergency services. It may be used for patients for whom supplemental oxygen is critical as long as the potential risks of introducing the device are assessed and additional risk control measures, where necessary, are applied. Risk controls may include monitoring and increased supervision of the patient.

There are a number of Firesafe™ Nozzle models available to cover a wide range of oxygen flowmeters, gas cylinders, oxygen concentrators and liquid supplies available in the market. Care should be taken to select the correct Firesafe™ Nozzle for the application.

The Firesafe™ Nozzle is intended for connection to reusable oxygen delivery devices, upstream of single patient use oxygen therapy sets and similar equipment. It is contraindicated for use downstream of aerosolised drug therapy devices. The device is non-sterile and has an intended life of 8 years.

## Warnings!

1. Read through this instruction before installing a Firesafe™ Nozzle. As with all medical equipment, attempting to use or install this device without a thorough understanding of its operation and limitations may result in patient or user injury.
2. This device offers a resistance to flow; consider the impact of this on the system when installing the device.
3. Do not store or install this device near an open flame or near a source of excessive heat that is likely to exceed that stated in the Device Specification.
4. Keep this device free from oil and/or grease. In particular, do not lubricate the hose barb connector.
5. Never administer oxygen or undertake oxygen therapy while smoking or when near an open flame.
6. Additional risk control measures (e.g. SpO2 monitoring with an alarm) are essential for patients who might suffer injury or death in the event that their oxygen flow is stopped unexpectedly.
7. This device is designed for use in the delivery of oxygen or oxygen enriched air. Do not use the device with any other gas.

## Cleaning, Maintenance & Disposal

Clean the exterior surfaces of the device using an alcohol or disinfectant wipe. Do not immerse the device in any fluid or allow fluid to enter the device. The Firesafe™ Nozzle is maintenance free and has an intended life of 8 years. Dispose of any Firesafe™ Nozzle 8 years after the manufacturing date engraved on the device. Once actuated, the Firesafe™ Nozzle cannot be reset and must be discarded.

## Device Specification

	≤ 1 kPa at 2 l/min
Resistance to Flow	≤ 2 kPa at 5 l/min
	≤ 16 kPa at 15 l/min
Flow Rate <sup>(1)</sup>	0.7 - 20 l/min
Maximum Operating Pressure <sup>(2)</sup>	1000 kPa
Maximum internal leakage following activation	10 ml/min
Maximum external leakage following activation	5 ml/min
Operating Temperature	0 to 50 °C (32 to 122 °F)
Transit & Storage Temperature range	-20 to 60 °C (-4 to 140 °F)
Humidity Range	0 to 100% RH

<sup>(1)</sup> The Firesafe™ Nozzle can be safely used at lower flow rates but may not meet the internal leakage rate specified when activated in the event of a fire.

<sup>(2)</sup> Maximum static pressure delivered by the gas supply source.

## Installation Instructions

1. Remove the existing connector from the oxygen supply source and fit the appropriate Firesafe™ Nozzle.
2. Check that gas flows freely through the Firesafe™ Nozzle and that there are no obvious leaks.

**Warning!** If a Firesafe™ Cannula Valve fitted downstream of a Firesafe Nozzle™ is activated by fire, the supply tube downstream of Firesafe™ Nozzle will rise until it reaches the maximum output pressure of the supply system. In order to ensure that the connection between the Firesafe™ Nozzle and the supply system can withstand this pressure, it is recommended that a type test be carried out with each style of connector or tubing used.

**Note:** When fitting to an oxygen concentrator, to meet the requirements of ISO 8359:1996 AMD1 (2012), the DISS version of the Firesafe™ Nozzle (REF 827-0032) should be tightened to a torque of 1.5 to 2.0 Nm (1.1 to 1.5 lbf.ft). This will ensure the Firesafe™ Nozzle can only be removed by the use of a tool. A specially adapted Firesafe™ Torque Wrench (REF 802-0017) is available for this purpose.

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BPR Medical Limited