

UDV01/AF – 13th June 2018

Re: Ultraflow Demand Valves – increase in intended life

Dear Customer,

In accordance with MHRA guidance 'Managing Medical Devices' (April 2015), the Medical Device Regulations and company policy, we are obliged to manage the risk of medical devices throughout their entire lifecycle. An essential principle of this is to establish a safe useful life for each device type and to communicate these to our customers.

There has been a significant change made to the lifespan and service interval of this device due to the consistent high quality and performance after a four year period.

We are pleased to announce that any demand valves ordered from 18th June 2018 onwards will have an extended lifespan of ten years instead of eight. This also means that the service interval is now five years instead of four.

The lifespan and service intervals for *serviceable* medical gas therapy products manufactured by BPR Medical are provided below for ease of reference.

<u>Product</u>	<u>Useful Life</u>	<u>Service Interval</u>
Dialflow Regulators (range A and E)	12 years	2 years
Dialflow Regulators (range C and D)	12 years	4 years
Dialflow Meters	12 years	4 years
Microdial Flowmeters	12 years	2 years
Pressure Regulators	12 years	4 years
OASIS	12 years	4 years
Ultraflow Demand Valves	10 years	5 years

If you are left with any questions regarding this notice, please contact me on alison.fisher@bprmedical.com or sales@bprmedical.com or call us on 01623 628 281 and we'll be happy to help.

Yours Faithfully



Alison Fisher
Sales Manager.